



ManuMix[®]

• Mischsystem für den klinischen Einsatz / Mixing system for clinical applications

ManuMix®

ManuMix® ist das patentierte sterile Mischsystem zum manuellen Anmischen von Knochenzement niedrig- bis mittelviskoser Konsistenz und nachfolgendem Transfer in Spritzen für den klinischen Einsatz.

ManuMix® is the patented sterile mixing system for bone cement of low to medium viscosity and subsequent transfer in syringes for clinical application.

- 1**
Mischen des Knochenzementes mit dem beigelegten Spatel
Mix bone cement with the included spatula*



- 2**
Fixieren der Spritze und Einsetzen des Siegels
Fix syringe and insert transfer device*



- 3**
Sauberer Transfer des Knochenzementes in Spritzen
Clean transfer of bone cement into syringes*



- 4**
- ▶ Verkürzte offene Produkthandhabung
 - ▶ Weniger Lufteinschlüsse
 - ▶ Sicherer und sauberer Transfer
 - ▶ Optimale Entnahmetechnik mit geringem Zementrückstand
 - ▶ Shorter open handling of product
 - ▶ Reduced air entrapment
 - ▶ Safe and clean transfer
 - ▶ Optimal removal technique with low cement residues



*Für Details bitte die Gebrauchsanweisung beachten
For details, please follow the instructions for use

ManuMix® QTY 10 Bestellnr./ Catalogue No.: 02-0602

responsible manufacturer:



aap Biomaterials GmbH
Lagerstraße 11-15 • 64807 Dieburg
Germany

Tel. + 49 6071 929-0
Fax + 49 6071 929-100
customer.service@aapbio.de
WWW.AAPBIO.DE